

# SEATING & TRANSPORT EQUIPMENT

## RIFTON MOBILE CHAIR



The High Back Mobile Chair, with push handle, pads, armrests and footboard.



Trunk Support adjustment is easily done by squeezing the blue clip at the back, slide along the bars then release.

Refer to page S - 23 for order information and options.



Push Handle & Locking Castors.



Flat or Contoured Headrest.



The Trunk Support strap snaps on and off.



Turning the handles under the seat can adjust the height or tilt easily.

### Low Back Mobile Chair Components/Options:

1. Tilting Chair Frame
2. Push Handle & lockable Castors
3. Angle Adjustable Low Backrest
4. Removable Abductor
5. Adjustable Footboard with Sandals

**NB:**

1. Must choose one frame.
2. Please choose backrest separately.



### High Back Mobile Chair Components/Options:

1. Tilting Chair Frame
2. Push Handle & lockable Castors
3. Adjustable Contoured Headrest
4. Angle Adjustable High Backrest
5. Trunk Support
6. Removable Abductor
7. Adjustable Footboard with Sandals



### Common Options:

- A. Tray provides working surface for different activities.
- B. Handhold helps child to maintain balance.
- C. Additional Side Pads can be added to reduce seat width and provide more comfortable support.
- D. Leg Prompts use the same attachment socket as an alternative to the Abductor. Straps closed just below the knee for improved security of legs.
- E. Arm Anchors as an alternative to the Arm Prompts for table activities. Available on right or left arm, in two sizes.

