


SEATING & TRANSPORT EQUIPMENT



The Advancement Chair can provide forward tilting, with footboard pushed back, encouraging sit-to-stand practice.



The Rifton Advancement Chair comes in three sizes. Push handle on small chair only.



Rifton Advancement Chair Accessories:

1. **Adjustable back and tilt frame** can provide functional or leisure support to the child and is able to adjust the degree of forward leaning for different activities and skill level of the user.
2. **Trunk Support (Small & Large)** can be adjusted in width, angle and height easily and can be lowered or removed if the child gains better trunk control.
3. **Arm Prompts** can be fitted onto any position of L-bar and is adjustable in height, rotating and tilting in space, giving the child upper body support and encourage head control.
4. **Swing-away L-bar** allows easy transfer of the child.
5. **Padded Groin Strap** assists in correct, comfortable pelvic positioning.
6. **Padded Abductor** is removable and depth adjustable.
7. **Sandals and Wedges** fix the feet on the footboard to improve midline position. The footboard slides out of the way for sit-to-stand practice.
8. **Locking Castors** allow chair to move around.
9. **Flat/Contoured Headrest** is height adjustable, offers comfortable support to the the head when it is needed.
10. **Leg Prompts** use same attachment socket and is an alternative to the Abductor. Straps closed just below the knee for improved security of legs.
11. **Arm Anchors** as an alternative to the Arm Prompts for table activities. Available for right or left arm, in two sizes.

