

LECKEY MYGO SEATING SYSTEM

Model: SOB-LKMG

"It's my go" is a common declaration of independence heard from children of ages 4 - 12. The Mygo Seating System from Leckey aims to give the child this freedom of mobility while upholding postural support in the pelvis, trunk and head, leg and foot positioning.

Features include:

- Unique leg support system that allows for shortening and wind sweeping up to 10° from the midline. The child's pelvis is encouraged to sit in a 'square' position.
- Flexible Sacral Support is a cushion that assists the child into a more extended posture. It spreads the support over a larger area, reducing pressure on individual areas by following the natural curve of the spine. The support is then adjusted to the position best for the child's postural support.
- Hip guides that work with the Sacral Support adjust in width so that the child will sit in the midline position each time.
- 4-point Positioning Harness encourages the pelvis to sit in a more horizontal alignment.
- The Upper Backrest Cushion is height adjustable. It has two lateral supports, adjustable in height, width and angle.
- The Cushioned Chest Support is a small padded cushion attached

to the lateral supports. The benefit of the small size means that the upper body and arm movement are not restricted whilst the trunk support is maintained.

- The Trunk Harness is a full vest-like cushioned support, giving postural assistance to the trunk and support to the broader chest area. It is secured into the Pelvic Positioning Harness.
- The Backrest Adjustment Feature of the Mygo Seating System allows the back of the chair to pivot at the anatomically correct position. With this feature, the depth and angle of the seat can be adjusted to accommodate for changes in posture, growth and seating angle over time.
- The Head Positioning works together with the Postural Support elements of the chair to further enhance and extend the support available for the child.
- Lower Leg Supports follow the Upper Leg Guides when moving laterally; however, when moving forward, the Lower Leg Supports will remain in the same position. The position of the child's feet, therefore, is always behind the knee line, thus giving better support during the child's growth period to improve posture.
- Height and Angle Adjustable Footplates: come with either sandal or ankle supports. Sandal supports allow for rotation, and angle-adjustment to cater for plantar-flexion and dorsiflexion.
- Cushioned Height and Angle Adjustable Armrests, easily removable for when elbow support is not needed.
- Cushioned Activity Tray.
- Mygo can fit onto the hi-lo indoor base or many other wheel chair bases.
- A child-friendly, colourful design.



Chest Support



Pelvic Harness



Mygo Seat on a Otto Bock Kimba Base.



Mygo Seat on a Hi-Lo Indoor Base.



Mygo Seat can be used to do all sorts of interesting activities.