

SEATING & TRANSPORT EQUIPMENT

E-Z-ON VEST

E-Z-ON MODIFIED VEST



The E-Z-ON Modified Vest allows a child to travel in a prone or supine position securely. This will suit certain disabilities such as a spinal injury or hip spica that requires the child to lie down during transport. The child can nap in the car. There is no need for modification/locking clips to be installed. Just put on the vest, lie the child on the seat, and secure the existing seat belts after threading them through the sides of the vest and the leg belt.




E-Z-ON Modified Vest helps to secure the child lying on the back seat of a vehicle.

Model No.	Age	Closure
SV101M2/XS	2 ~ 5 yrs	Double D Ring
SV101M2	5 ~ 12 yrs	D Ring/Optional Push Button



E-Z-ON Modified Vest comprises a vest with crotch straps and a leg belt.



E-Z-ON Modified Vest allows transportation of a child in a prone or supine position.

E-Z-ON PARTS & ACCESSORIES

- A. Sherpa Neck Pad - SV100NP
- B. Crotch Straps - SV100CRP
To retrofit a vest that does not come with crotch strap.
- C. Zipper Inserts - SV100AZ
To adj. vest size by expanding zipper. 4 sizes available for XXS, S, M, & L vest.
- D. Wrist Restraints - V100WRP
- E. Ankle Restraints - SV100AR
- F. Locking Clip - SV4012
- G. Webbing Cutter - SV1004
- H. Upper Arm Restraints - SV100UA

