

RIFTON PRONE STANDER WITH TRAY

The Rifton Prone Stander is easily and accurately tilted with a hand crank from horizontal to vertical (85°). Removable lateral support blocks can be moved to any position along the length of the main board and at an angle up to 30°. Adjust the main board and the footboard to the height best suited to the user. The knee board in AE930 is removable for use by smaller children and infants. When the Stander is in a vertical position the footboard is at almost floor level for easier transfer. You can also tilt the stander to horizontal for easier transfer and set up. Non-marking 75 cm double-locking castors prevent roll and swivel. Optional knee lateral pads can be adjusted vertically, horizontally and rotated to provided better knee positioning. All models include a tray.

Optional accessories include Abduction Block, Abduction Wedge and Hip Stabilizer, Suction Cup Handholds, Sandals and Wedges, etc.

Model No	AE930	AE940	AE950
User Height, cm	64 - 122	112- 150	145 - 182
Max. working load, kg	45	68	91
Base W x L, cm	61 x 89	66 x 99	76 x 109
Length of board, cm	41 or 60	76	99
Distance between Lat. Supp., cm	20 - 41	20 - 41	30 - 56
Top of board to footboard, cm	46 - 89	84 - 110	107 - 140
Board width (top/bottom), cm	32/32	28/41	38/48
Tray (included) size (WxL), cm	57 x 44	57 x 44	75 x 44
Round Abduction Block, 10cm D.	AE933	AE933	AE933
Rd. Abduc. Block w/hip Stabilizer	AE655	AE655	-
Collar for AE933, add 5 cm	AE934	AE934	AE934
Oblong Abduction Block	AE653	-	-
Pr of Suction Cup Handholds	AK821	AK821	AK821
Pr of Sandals	<-- AE422 or	AE433 -->	AE433 or AE432
Sandal size, cm	(S-8x18)	(M-10x23)	(L 13x27)
Pr of Foot Wedges (Size)	<- AE424 (S) or	AE435(M) -->	AE435 or AE434(L)
Adj. Abduction Wedge	-	AE945	AE945
Adj. Abd. Wedge w/hip Stabilizer	-	AE605	AE605
Extra pr of Trunk Lateral Support	-	AE614	AE614
Pr of Knee Lateral Support	AE942	AE942	AE942



Prone Stander tray provides enough space for various activities.



Footboard can adjust to suit smaller children.



Round Abductor w/hip stabilizer keeps leg apart and improves pelvic position.



Knee Lateral Supports control knee position.



Hand crank winding with position indicator makes it easy for the caregiver to position the user comfortably.



Sandals & wedges aid foot positioning and body alignment.