

# Adaptive Equipment

## T-Tables

Often a therapist requires a surface to facilitate upper extremity movements once a child has been correctly seated. The T-tables are ideal for this purpose as they are lightweight and portable. The T-tables have been designed to provide a tilting surface adjustable between 0° (horizontal surface) and 90° (fully upright surface) to support toys or educational materials. They can be used with various equipment including Benches, Kinder Chairs, *Tumble Forms* Floor Sitters, Foam Logs, Balls, Standing Frames and wheelchairs.

Two adjustable clamps can hold books and paper. The wooden ledge can be used to prevent items from sliding off the table, or removed to provide a comfortable forearm position.

Model No	<i>AT1</i>	<i>AT2</i>
Size	Small	Large
Table surface size, cm	37.5 x 51.5	50.5 x 75.5
Height adj. cm	30 - 43	45 - 64.5
Width under table top	37	55
Surface tilts (degree)	0 - 90	0 - 90



T-Table AT2 works well with a therapy bench (see page S-13) to provide a play surface for the child. These tables help prepares the child's shoulders, arms and hands for self help skills and school readiness.



T-Table AT1 and AT2 can provide surface tilting from 0 to 90 degrees to suit different activities.

## Dynamic Standing Frame

This standing frame is designed to provide balance and allow limited movement in all directions. It features an adjustable heel bar, a vertical adjustable knee support, a thoracic or pelvic padded strap and a non-slip rubber mat. Optional extra padded lumbar restraint, is available for the Medium and Large Standers.

Model No	<i>ADF9314</i>	<i>ADF9315</i>	<i>ADF9316</i>
Size	Small	Medium	Large
Height adj., cm	45 - 60	65 - 85	80 - 105
Base Dimension, cm	43 x 60	43 x 60	59 x 79
Options:			
<i>Extra Lumbar Restraint</i>	<i>ADF93141</i>	<i>ADF93141</i>	<i>ADF93161</i>



# Please note that the Activity Tray has been discontinued

