

## LECKEY HORIZON STANDING SYSTEM

The Leckey Horizon Stander is a power assisted standing frame helping to assist with many of the current lifting and handling issues that therapists and caregiver frequently experience. The standing frame comes with a impressive level of postural support giving the user a high level of comfort while using the frame. Many of the supports are adjustable to contour to the users body shape.

The Leckey Horizon comes in two sizes and the frame can be used in either a **Supine** or **Prone** position. Converting the stander from a Supine to a Prone position can be done in three easy steps. Simply remove the adjustable headrest, replace the supine tray with a prone tray and reverse the sandals. It is that easy.

### Standard Features:

**Head Support** is multi adjustable (Supine Stander only).

**Padded Back/Chest Support** adjusts in height to allow for growth.

**Lateral Supports** are flexible and adjust in height and width.

**Activity Tray** is adjustable.

**Hip Support with Hip Belt** is height depth and angle adjustable.

**Knee Cups** are depth and angle adjustment.

**Swivel Castors** are lockable.

**Footplates** are Plantarflexion and Dorsiflexion.

**Sandals** are multi adjustable.

**Acuator and Headset.**



Leckey Horizon Standing Frame in Prone position with an angle adjustable tray.



Leckey Horizon Supine Standing Frame.

**Supine Stander (Conversion) Kit** (AOBLSK1 - Size 1)  
(AOBLSK2 - Size 2)

**Prone Stander (Conversion) Kit** (AOBLPK1 - Size 1)  
(AOBLPK2 - Size 2)

### Optional Accessories

**Lateral Head Supports**, height and width adjusts for Supine Stander.

**Lateral Hip Supports**

**Whitmyer Pro Series Head Support**

## MINOR SUPINE STANDER

**Model No. :** **AMSS-E Electric Tilt**  
**AMSS-G Gas Tilt**

The Minor frame is made from light weight aluminium alloy. This standing frame helps develop a child's independence, dignity and self esteem while at the same time promoting muscle and bone development for children with conditions like cerebral palsy, paraplegia and hemiplegia.

### Main Features:

- Components adjust along the centre column to suit individual user.
- Horizontal top for easier transfer and positioning.
- Folds for easy storage.
- Electric or Gas operated tilt system.
- Four lockable swivel castors.
- Removable therapy table.
- Single footrests

**NEW PRODUCT**



Minor Supine Stander

Leckey Horizon Stander	AOBLH1	AOBLH2
Size	1	2
Age (approx.)	4 - 10	10 - 18
Length, cm	113	128
Width, cm	58	64
Horizontal Height , cm	78	80
Product Weight, kg	41	47
Maximum User Weight, kg	50	80
User Standing Height (Min to Max), cm	100 - 140	140-180
Shoulder Height (Min to Max), cm	82 - 115	100 - 145
Hip Pad Height (Min to Max), cm	45 - 75	75 - 96
Knee Height (Min to Max), cm	25 - 45	40 - 55
Width between Chest Laterals, cm	20 - 30	24 - 38
Footplate Angle Adjustment	10°	10°
Tray Sizes (L x W), cm	42 x 44	52 x 59
Tray Angle Adjustment (in Prone)	0° - 40°	0° - 40°
Tray Angle Adjustment (in Supine)	20° - 30°	20° - 30°