


# BATHING AIDS



Blue Wave Bath Chair works with or without stand.

The Blue Wave Bath Tub Stand will raise the seat height to 43 cm.

The Blue Wave Shower Base raises the seat to a height of 71 cm, allowing the user to be wheeled from one room straight into the shower recess.



Bath chair can be used in lounge room, garden or beach.



Standard fabric (vinyl-covered nylon) provides a firm surface for postural support.



Height can be adjusted with child in chair.



Soft fabric (polyester knit) for those who prefer sink-in comfort rather than postural support.



Seat and back angle adjusts easily by triggers.



A smaller bath chair can be converted to larger size with conversion kit.



Fabric can be removed for cleaning.

Size	Small	Medium	Large
<b>Blue Wave Bath Chair (standard fibre)</b>	<b>BE514</b>	<b>BE524</b>	<b>BE534</b>
<b>Blue Wave Bath Chair (new soft fibre)</b>	<b>BE513</b>	<b>BE523</b>	<b>BE533</b>
<b>Blue Wave Tub Stand</b>	<-----	<b>BE517</b>	----->
<b>Blue Wave Shower Base</b>	<-----	<b>BE518</b>	----->
Height of user, cm	up to 91	81 - 127	117 - 173
Max. working load, kg	27	54	73
Back height, cm	46	64	81
Seat Depth x Width, cm	33 x 37	33 x 37	42 x 37
Seat Height, cm	5 - 18	5 - 18	5 - 18
Seat ht. on Tub Stand, cm	<-----	43	----->
Tub Stand W x L x H, cm	<-----	43 x 81 x 25	----->
Seat ht. on Shower Stand, cm	<-----	71	----->
Shower Stand size WxLxH, cm	<-----	69 x 84 x 53	----->
<b>Conversion Kit, BE514 to BE524</b>	<b>BE519</b>	<i>n/a</i>	<i>n/a</i>
<b>Conversion Kit, BE524 to BE534</b>	<i>n/a</i>	<b>BE527</b>	<i>n/a</i>