



(FAS Therapeutic Equipment P/L) Wonderland Rehab & Child Care Products

Newsletter July, 2008

Call Us For Free Trial & Quote!

In this month newsletter, we would like to introduce Kidster Wheelchair and Swivel Base for Columbia Car Seat in detail. You may find them useful for your clients.

Kidster Wheelchair

Kidster Wheelchair is a grow-with-child manual tilt-in-space wheelchair, it comes with two sizes, 10-13 inches and 14-20 inches. The frame base is adjustable for growth in term of seat width, seat depth and back height. The base allows to turn around for rear wheel/front wheel driven. Dynamic Recline is from 90-180 degree. Posterior tilt can go up to 80 degree. Full customized seat and back cushion will be suitable for all different condition clients. The latest concept can be integrated in this chair as well, lateral tilt, which allows the seat to tilt horizontally for providing more options to the seating posture. More options and frame colours are available.



Rear Wheel Driven



Tilt in space
(up to 80 degree)



Customized seating system with options
- Multi-point head rest
- Swing away lateral support
- Hip pads
- Pommel
- Harness and more



Front Wheel Driven



Dynamic Recline
(from 90-180 degree)



Growable frame



Anti-tipper & Transport Anchor



Lateral Tilt

Swivel Base for Columbia Car Seat

Have you found a lot of trouble while transferring a child from car seat to stroller or vice versa?

Now, here is a solution, Swivel base for Columbia Car Seat, just pull couple pins on the base, it allows you to turn the car seat in 90 degree, facing car door, for easier transfer. Only available on Columbia Car Seat (Child size).



Wonderland Rehab & Child Care Products

Tel.: 02 9649 4044 (NSW)
1300 30 35 36 (local call)
Fax: 02 9649 4055 (NSW)
Email: wonderland@adsl.on.net

FAS Therapeutic Equipment P/L

Tel. : 03 9587 6766 (Other States)
1300 30 35 36 (local call)
Fax: 03 9587 6899 (Other States)
Email: sales@fasequipment.com

Website: www.fasequipment.com